97050173
rev. 014
1 7 / 0 5

CE
0051

 myray

WDS X-POD

IT
EN
FR
DE
ES
PT
PL
RU
ZH
CS
SV

INDICE

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A. FCC Regulatory information

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This device contains FCC ID: POOWML-C40.

CONTENTS

1 Introduction	2	13 Possible operating modes	18
1.1 Conventions in text.....	2	13.1 Stand alone mode.....	18
2 Safety guidelines	2	13.2 Continuous connection mode.....	18
2.1 Warnings for use.....	3	14 Handheld WDS X-POD	19
2.2 Safety and Health Care.....	3	14.1 Fast reference guide.....	19
2.3 Maintenance and Disposal.....	4	14.2 Firmware update from USB®.....	19
2.4 Returning the equipment.....	4	14.3 Detailed handheld operation.....	19
3 Warranty	4	14.3.1 Setting up the handheld.....	19
3.1 Software not covered by warranty.....	5	14.3.2 Introduction.....	20
3.2 Battery warranty limitations.....	5	14.3.2.1 Charging and discharging.....	20
3.3 Limitation of liability.....	5	14.3.2.2 Turning the sensor on and off.....	22
4 Standards	5	14.3.2.3 SD Card.....	23
5 Box content	6	14.3.3 Capturing images.....	24
6 Product Identification	8	14.3.4 Managing and consulting the images.....	26
6.1 Intraoral sensor.....	8	14.3.4.1 Editing and measuring images.....	28
6.2 Handheld WDS X-POD.....	8	14.3.4.2 Modifying image orientation.....	30
6.3 Battery charger (power supply).....	9	14.3.5 Settings menu.....	31
7 WDS X-POD system features	10	14.3.6 Indicator LED.....	33
7.1 Environmental characteristics.....	10	14.3.7 Image transfer.....	33
7.2 Handheld characteristics.....	10	14.3.8 Handheld display screen.....	34
7.3 Sensor characteristics.....	11	15 Quality of the x-ray images	34
7.4 Battery charger specifications:.....	11	16 Problem solving	35
7.5 Battery characteristics.....	11		
7.6 Handheld mount.....	11		
7.6.1 Securing the support to the x-ray unit arm (post).....	12		
7.6.2 Installing the wall mount.....	12		
7.7 Personal computer.....	13		
7.8 Bluetooth.....	13		
7.9 SD Card.....	13		
7.10 USB®.....	14		
8 Compatibility with x-ray generators	14		
9 Disposable covers	15		
10 Cleaning and disinfecting	15		
11 Material preparation	17		
12 Use of the centering device	17		

1 Introduction

WDS X-POD is a MyRay® portable medical device capable of immediately capturing diagnostic quality intraoral x-ray images and display them on the handheld's high resolution color LCD screen.

Along with viewing the x-rays, it is possible to perform some simple tasks like zooming, adjusting brightness, range, contrast, etc.. by tapping the touch screen with your bare fingers or when wearing gloves.

In addition, it is possible to:

- transfer the images to a Personal Computer (PC), while they are being captured, via Bluetooth® wireless protocol
- collect hundreds of images on the SD memory card (Secure Digital), organizing them in named folders or identifying them with descriptions.
- transfer the images to a PC with a fast USB® connection or SD card reader connected to the PC at a later time.

After the images have been transferred to a PC, they can be viewed on a large screen and permanently saved in the patient folders using any dental image management application.

1.1 Conventions in text

The following abbreviations will be used:

PC = Personal Computer

X-POD = MyRay system WDS X-POD

CEFLA sc = CEFLA sc, Imola (BO), Italy – also referred to as the MANUFACTURER

RF = Radio frequency

LAN = Local Area Network

LED = indicator light (Light Emitting Diode)

MyRay = registered trademark of CEFLA sc, Imola (BO), Italy

Microsoft®, Windows®, Apple®, Bluetooth®, USB® are registered trademarks. They will be referred to in the text only in reference to the technology and work method that these trademarks represent.

2 Safety guidelines

Pay special attention to the sections of the manual in which the indicated symbols appear:

	Warnings regarding operator or patient safety
	Warnings regarding risk of product damage or failure to operate, important warnings regarding the warranty.
	The manufacturer's website contains a list of authorised agents.

- The X-POD system and relative software drivers are developed and manufactured by CEFLA sc.
- The instructions inform the user on how to properly operate the X-POD system. Carefully read and become familiar with the content of this manual before attempting to use the equipment.
- This product is intended only to capture and display intraoral x-ray images therefore the diagnostic study is left entirely up to the dental professional.
- In order to use the X-POD system with a PC a software application for managing and storing the dental x-rays is required.
- Consult the relative manual supplied by the software manufacturer for instructions on how to install and use the image management software.
- It is forbidden to reproduce, store and transmit this publication by any means (electronically, mechanically, by photocopying, translating or in other ways) without a written authorisation issued by the Manufacturer.
- CEFLA sc maintains a company policy of constant product development and improvement therefore some specific instructions and images found in this manual may differ from the product actually purchased. CEFLA sc reserves the right to make any changes to this manual without giving prior notice.
- The information, technical specifications and illustrations contained in this publication are not binding. CEFLA sc reserves the right to make technical improvements and changes without modifying these instructions.
- Important: in accordance with privacy laws in force in several countries, all sensitive personal information must be adequately protected. In addition, patients must sign a consent form before personal information or images are transmitted across networks.
- The original version of this manual is written in Italian.0
- Translation from the original in Italian
- All the registered trade marks and product names mentioned in this manual are the property of their respective owners.

2.1 Warnings for use

- Electronic apparatuses may cause or be subject to interferences when used near other electromagnetic equipment such as mobile phones, personal computers equipped with Wireless LAN cards and microwave ovens. Keep the parts of the ZX-POD system and Personal Computer used to acquire and save the x-ray images away from RF sources such as wireless LAN cards, other radiation devices, Home RF devices, microwave ovens; the recommended distance is at least 2 meters in the case of microwave ovens and 1 meter in all other cases.
- We recommend to regularly make backup copies of all acquired images.
- If the personal computer (optional) which the images are transferred to is connected to Internet an updated antivirus program is highly recommended and the PC should be used only for work purposes.
- **Important information:** if the PC fails while the x-ray images are being transferred (software or operating system crash, personal computer failure) the images are still saved in system memory X-POD. An SD memory card is used, therefore the images are not lost and can either be transferred at a later time or recovered by extracting memory data and reading it with an SD card reader.
- Consult the manuals that deal with the individual devices (PC, x-ray unit, printer, etc.) to use the equipment together with the X-POD system.
- Only specially trained technicians should install the other system components (computer or computer network, software to manage and save images, x-ray generator, etc...). In particular, keep in mind that installation of the x-ray equipment is checked and inspected by a qualified technician.

	<p>The USB connection[®] used by the X-POD system is not a normal electrical connection but requires specific cables for USB[®] 2.0 (marked USB[®] 2.0 Hi-speed[®], see paragraph "Product identification"). the cable is type "A → micro B" (supplied).</p>
---	--

	<p>The sensor used to capture images is fragile and bothered by static electricity. Handle it with care. Do not deform or squeeze it.</p> <p>Do not touch electric contacts when the connector is not inserted in place on the interface.</p> <p>Do not disconnect the sensor while it is turned on without first briefly pressing the power on/off button to shut off electrical power to it. As a result, the sensor can be safely removed even if careful attention is not paid. See paragraph "Using the sensor". When reconnected, power will automatically be supplied to the sensor.</p>
---	---

2.2 Safety and Health Care

- X-POD is a medical device for intraoral x-ray acquisition. It is intended for use only by qualified dental professionals. Do not use the system for any purposes other than capturing intraoral radiographs. The device can be used only by clinicians well instructed in the fields of dentistry and radiology.
- Do NOT use X-POD in the vicinity of life support (e.g.: pacemaker or cardiac stimulators) and hearing aids even if the electromagnetic field emitted by the sensor is unremarkable. In medical facilities, before using any electronic devices, always assure they are compatible with the equipment present.
- **The disposable protection control sheaths provided are mandatory to safeguard the patients against cross-contamination. The disposable protection control sheaths are class IIa medical devices. They cannot be replaced with others that are poorer quality. Contact the dealer who supplied the sensor to obtain additional disposable protection control sheaths.**
- Cover all the components which personnel may touch with disposable infection control sheaths as they may be contaminated through indirect contact with the patient's mouth. In particular, be careful when handling the PC mouse, keyboard or touch screen and X-POD.
- Some parts (USB cable[®], silicon rubber protection, disposable infection control sheath, items used for packaging, x-ray sensor) may result in suffocation if swallowed or used erroneously and inappropriately. Avoid unintended, inappropriate and misuse and keep out of reach of children.
- Use only the battery charger provided that conforms to I.E.C. 60601-1 IIa ed.

- When the images are synchronized with the USB port[®], NO x-rays can be taken. Do not leave the sensor inside the patient's mouth.
- The handheld is to be connected to a PC with a USB cable[®]. In this case, the PC should comply to IEC 60950 and its connection to the power mains must have an additional ground connection or isolation transformer. The system (device + PC) complies to IEC 60601-1-1 in this configuration.
- Never use in the presence of mixtures of flammable anesthetic gas with air, oxygen or nitrogen protoxide.
- Installation must meet the essential requirements of I.E.C. 60601-1-1 (General requirements for safety).
- The electric installation in the room in which the apparatus is installed must comply to standard IEC 60364-7-710 (Standards regarding electrical installations in medical locations).

2.3 Maintenance and Disposal

The equipment does not contain any parts that can be repaired directly by the user. The internal battery needs to be replaced by a specialized technician at a Service center when it is dead. In the event of malfunctions, do not try to perform any maintenance jobs. Contact the manufacturer or local dealer directly at the phone numbers given on the warranty certificate. If the apparatus has to be returned to the manufacturer or Service Center for any reason, completely disinfect the outside of the apparatus with a specific product (see paragraph "Cleaning and disinfecting") and send it back preferably in its original box.

No electronic parts of the X-POD require maintenance. Devices may be broken and the warranty is void if the casings of the handheld, sensor, battery charger, SD card and Bluetooth dongle[®] are opened to reach the circuits inside.

Do not use the sensor on a patient if a system malfunction is present or suspected.

Discard the apparatus in accordance with current regulations at the end of its service life.

	<p>Always disinfect all the parts outside the equipment before disposing of it. Sort the materials for separate disposal.</p>
--	---

In compliance with Directives 2011/65/EU - 2012/19/UE regarding restriction of the use of certain hazardous substances in electrical and electronic equipment along with waste electrical and electronic equipment, it is forbidden to dispose of this equipment in the municipal waste stream as unsorted municipal waste. This equipment must be sorted and disposed of separately. When new equipment that is similar is purchased, the old equipment must be given to the dealer for disposal. As regards reuse, recycling and other forms of recovery of waste electrical and electronic equipment, the producer carries out the functions set forth by laws and regulations currently in force in the individual countries. A high level of separate collection of waste electrical and electronic equipment is indispensable to efficiently recycle, treat and dispose of the equipment. Recycling and treatment operations should comply with minimum standards to assure human health and high environmental protection as well as favor recycling of the materials included in the equipment. The symbol indicating separate collection for electrical and electronic equipment consists of the crossed out bin marked on the equipment. Illegal waste clearance and disposal shall be punished as established by laws and regulations currently in force in the individual countries.

2.4 Returning the equipment

If the equipment needs to be returned to the manufacturer or Service Center for any reason whatsoever, follow the disinfecting and safety instructions provided (see paragraphs "Cleaning and disinfecting" and "Maintenance and disposal") and do not send the SD card in order to assure security and privacy of the data and images stored on it.

In order to protect the privacy rights of patients and operators in compliance with privacy laws, if an SD card is received, the manufacturer shall not be held responsible for its content and will format it as a preventive measure.

It will not at all be possible to recover any images stored on it.

3 Warranty

The manufacturer warrants safety, reliability and performance of the equipment. The warranty is valid provided that the following terms and conditions are observed:

- Observing the conditions given on the warranty certificate.

The appliance must be used exclusively in accordance with the instructions contained in this instruction manual.

The electric installation in the room in which the apparatus is installed must comply to standard IEC 60364-7-710 (Standards regarding electrical installations in medical locations).

Installation, repairs, settings and in general all the operations for which the equipment covers need to be opened must be carried out only by technicians authorized by the manufacturer.

WARNING!

The manufacturer shall not be held liable for any personal injury or property damage arising from failure to heed the following clauses.

3.1 Software not covered by warranty

The Software is supplied in its original condition and the Manufacturer shall not be held liable or warrant any original defects or defects originating over time and shall not guarantee quality and proper operation of the software. In addition, the manufacturer shall not honor or provide any warranty regarding conformity of the software to the information given on-line or in electronic documentation or in any case made available except for the warranty for the CD-ROM if damaged or unusable.

Any warranty is also excluded for Software integrated in - or otherwise being a part of - other Software applications developed by third parties. As far as these applications are concerned, the Manufacturer also expressly declares not to have carried out and not to carry out any inspection activity or other activities to guarantee the software operation.

3.2 Battery warranty limitations

The battery provided in the X-POD system is covered by a 6-month warranty in effect as of the date of installation.

3.3 Limitation of liability

Under no circumstances shall the manufacturer or its suppliers be held liable for any direct or indirect damages (including damages for no or lost earnings or savings, interruption of practice, loss of information or data and any other economic losses) faced by the users, i.e. third parties, deriving from use or failure to use the software, even in the event that the manufacturer was informed of the possibility of said damaged.

Liability is limited not only if the software is not used as instructed by the manufacturer but also in the event of unintended and erroneous use.

4 Standards

X-POD is a class IIa medical device for taking intraoral x-rays.

X-POD and its accessories are designed and constructed in compliance with the following standards:

92/31/CEE	European electromagnetic compatibility directive.
93/42/CEE and subsequent changes	Apparatus classified as class IIa in accordance with the Medical electrical equipment directive.
EN 60601-1:1998 - IIa ed.	Safety standards for medical electrical equipment. Apparatus designed to work under environmental conditions present in closed work areas in accordance with the parameters set forth by these standards.
EN 60601-1-2:2003	Medical electrical equipment: Electromagnetic compatibility - Requirements and tests.
EN 60878:2003	Graphic symbols for medical electrical equipment functions.
UNI CEI EN ISO 14971:2007	Medical devices. Medical device risk management application.
1999/5/CE	Radiation equipment Directive LV standard.
ETSI EN 300 328-2 V1.6.1 2004	Radiation-emitting device regulation.
ETSI EN 301 489-1 V1.4.12006	Radiation equipment: collateral standard for electromagnetic compatibility.
ETSI EN 301 489-17 V1.2.12002	Radiation equipment: collateral standard for electromagnetic compatibility.
FCC part 15 subpart C	USA regulation for radiation-emitting devices not subject to end user license agreement.

The device contains a Bluetooth class I transmitter (Module FCC ID: POOWML-C40).
 Class II apparatus with type B applied part in accordance with standard I.E.C. 60601-1.
 Not protected against the penetration of liquids (degree of protection (IPX0)).

Never use in the presence of mixtures of flammable anesthetic gas with air, oxygen or nitrogen protoxide.

Do not install in areas exposed to freezing temperatures.

The apparatus complies with European and North American Standards for safety of electrical equipment and electromagnetic interferences therefore dangerous electromagnetic interferences are not generated and the apparatus is not affected by emissions released during operation of other equipment.



5 Box content

The package contains three tiers. Be careful not to discard the package if there are still components inside it.


Package components:

	<p>Handheld WDS X-POD</p>
	<p>MyRay sensor[®] with cable and connector</p>
	<p>MyRay centering device introductory kit, including disposable infection control sheaths.</p>
	<p>CD-Rom with the software and drivers</p>
	<p>Battery charger</p>
	<p>USB[®]2 HighSpeed[®] extension</p>
	<p>Bluetooth receiver[®] (Dongle USB[®])</p>
	<p>Hub USB</p>

	Warranty certificate
	SD Card
	Compliance certificate
	Belt
	Mount
	Neck piece
	Kit of power supply plugs

6 Product Identification

6.1 Intraoral sensor

The manufacturer's brand name (logo) and serial number of the part are given on the back of the sensor.

	MyRay logo [®] , registered trademark of CEFLA sc, Imola, Italy
W1234-123	Product identification and serial number of part: - Wxxx-xxx MyRay sensor size 1 - Wyyy-yyy MyRay sensor size 2

6.2 Handheld WDS X-POD

The brand name (logo) of the company that markets the sensor is given on the front. Further information is given on the warranty certificate.

	MyRay [®] trademark
	Name of the appliance

The nameplates (laser marking) are attached to the edges of handheld WDS X-POD.
The following information is printed on the plastic:


CEFLA sc	Manufacturer's name
WDS X-POD	Name of the appliance
WX00500	Serial number
1.5 A	Rated current and voltage

	ETL trademark (Control Number 3194988)
FCC ID: POOWML-C40	Device identification in accordance with FCC (it may not be present in devices not intended for USA / Canada markets)
	"Type B part attached" symbol in accordance with standard EN 60601-1
	Important information, consult the enclosed documentation
CE 0051	CE mark of approval
	Symbol IEC5010 for power on/off button
	"Read the instructions manual" symbol
	"Important information" symbol
	Ukraine compliance mark

6.3 Battery charger (power supply)

	<p>Use only the FRIWO brand model MPP 15 power supply provided code FW 7555M/09. Use of any other power supply may reduce safety of the device and put the patient or operator at risk.</p> <p>IMPORTANT INFORMATION: The battery charger is approved for use in North America (USA) or Canada and can run only on 120 VAC.</p>
--	--

The following information is stamped on the battery charger. The following information is given on the nameplate:

<p>Wall mount power supply labelling 64300006:</p> 	<p>(hot stamped on the plastic canister side)</p>	Manufacturer Logo
		Manufacturing Country
		Product name
		Plate data
		Symbol "Caution - Read the instruction manual"
		WARNING
		European Markings
		Rated current and voltage
		Production batch number. The production batch is given in format MMY, for example 0409 means the unit was produced in April 2009.

7 WDS X-POD system features

The technical features and some of the main system functions largely depend on the characteristics of the x-ray generator and software used to display and store the images.

A constant potential x-ray generator (DC) with long rectangular collimator (skin - focal spot distance not less than 30cm) should be used to obtain the best possible results.

Old generation x-ray units that do not allow the exposure times to be sufficiently reduced may not be suitable for use with the X-POD .

The equipment is designed to work under indoor environmental conditions in compliance with the parameters set forth by standards CEI EN 60601-1.T.

7.1 Environmental characteristics

Operation :

- temperature +10°C - +40°C.
- relative humidity ranging from 20% to 90%
- atmospheric pressure 500 - 1060 hPa

Storage:

- temperature -10 °C - +70 °C
- relative humidity ranging from 0% to 100%
- atmospheric pressure 500 - 1060 hPa

Battery storage:

- temperature -20°C - +45°C (up to 1 month) / -20°C - +35°C (up to 6 months)
- the apparatus is suitable for continuous operation

Neither the sensor (handheld), sensor connector or battery charger are protected against the penetration of liquids (IPX0 degree of protection).

7.2 Handheld characteristics

- Power supply: with rechargeable battery, Li-Poly, 3.7V, 5200 mAh
- PC power supply: 5V DC 500mA max, supplied by personal computer USB port ®.
- Fast charger: 9V DC 1500mA max
- Dimensions: 142 x 83 x 31mm (LxHxT), including connector area 154 x 99 x 31mm
- Weight: 380g
- Display/Touchscreen: 4.3" (95x54 mm); 480x272 pixel; 16.7 million colors/ Anti-glare; white led backlighting, 500cd/m²
- SD-Card: SD / SDHC compatible; 2 Gb SD-Card supplied
- USB®: USB® 2.0 HighSpeed / full speed
- Wireless: Bluetooth 2.0 EDR, class I
- Processor: 533 Mhz, 32 bit
- RAM: 64 Mb DDR
- Native image: acquisition at 14 bit – max. 16384 grey levels
- Image format for export: 16384 grey levels, Microsoft® Windows compatible

7.3 Sensor characteristics

MyRay sensors have a Cesium Iodide scintillator with layer of fiber optics

The front corners of the sensitive area of the sensor are smoothed and the protective shell has rounded corners for improved ergonomics and patient comfort.

- Pixel dimensions 20x20µm;
- Theoretical maximum resolution 25 lp/mm
- High sensitivity and dynamics: saturation 40mR normal, 14 bit converter 58 db dynamic range
- CMOS silicon layer life: more than 100000 exposures at 60 KV/4 mSv.
- CMOS protected against direct x-ray exposure with a layer of fiber optics.
- Sensor cable length: 850/1350 mm

Two intraoral sensors are available:

- MyRay sensor size 1:
 - External dimensions: 38.9 x 24.9 x 5.3 mm (Length x Width x Thickness at edges)
 - Active area: 30 x 20 mm
 - Pixel matrix: 1500 x 1000
- MyRay sensor size 2:
 - External dimensions: 41.9 x 30.4 x 5.7 mm (Length x Width x Thickness at edges)
 - Active area: 34 x 26 mm
 - Pixel matrix: 1700 x 1300

The cable is connected on the back of the sensor with a 13 mm thick coupler.

7.4 Battery charger specifications:

- Medical device battery charger 9V - 1.5 A model FRIWO MPP 15 code FW 7555M/09
- Mains supply voltage: 100 – 240 VAC, 50/60 Hz, 400 mA
- For USA and Canada: 120 VAC nominal
- Kit of universal sockets (EU, UK, USA, AUS, ROW)
- Power plug diam. 4 mm, central positive

7.5 Battery characteristics

Lithium-polymer batteries are installed inside the X-POD. Different kinds of similar batteries are used:

- Accu Italia, model 02820002 LCP 10/52/110
- Guangzhou Markyn Battery Co. Model 9051109

Only properly trained personnel are authorized to change the batteries.

7.6 Handheld mount


It comes with a wall or post mount useful to mount the WDS X-POD handheld on a cabinet or arm of the x-ray unit being used. It is also possible to install on commercial mounts RAM-A and RAM-B available from the Ram Mounts catalogue (www.ram-mount.com).

- Mounting plate: 84 x 59 mm (post or wall mounting)
- Mounting: silicon band for installation on post. Set up for expansion plugs for wall mounting
- Handheld mounting: fast, 360° rotation
- The handheld can be turned 360° without taking it off its mount; it can be removed by raising it with your hand.

The support bar needs to be secured if the handheld is to be used together with the support.

- Attach the support bar to the handheld making sure the direction is correct.


Place this part on the side where the power key is found.

7.6.1 Securing the support to the x-ray unit arm (post)

Perform these operations to secure the support to the x-ray unit arm, or more generally to a post:

- Insert the silicone band on the handheld support (see figure below) and secure around the desired post.
- Attach the safety lock to fasten the silicone band.


Lastly, place the device in its compartment.

7.6.2 Installing the wall mount

Perform these operations to install the wall mount:

- Remove the silicone band, if it has already been placed on the handheld support, and take off the two plugs (2).
- Drill a hole in the wall (8 mm diameter) at the height at which you want to hang the handheld. When drilling, be careful not to puncture the tubes or electrical cables found inside the wall.

Pos.	Code	Denomination	Qty
1		MYRAY TIE SUPPORT	1
2	97460387	CYL. PLUG D 13	2
3	56004017	PLAST. INSERT D 8 TPF4/V WITH SCREW ZN YELLOW /HOLE IN WALL D 8mm	1

- Move 25 mm down maintaining alignment with the first hole and make another hole that has the same diameter.
- Place the two inserts (3) supplied in the holes.
- Attach the handheld support (1) and tighten the two screws provided.
- Put the white plugs (2) back in place.
- Lastly, place the handheld inside the compartment provided.

7.7 Personal computer

(Not supplied)

Operating system:

- Microsoft® Windows® XP - Service Pack 2 or 3
- Microsoft® Windows® Vista
- Apple® Mac OS X 10.5 Leopard

Processor and RAM: according to the minimum operating system requirements.

Space on hard disk: 150 MB minimum.

Video card resolution: 1024x768 minimum with 16 million colors (32bit true color).

Other peripherals:

- CD-ROM
- Integrated bluetooth or with a USB adapter (dongle)

7.8 Bluetooth

X-POD uses a class I Bluetooth communication channel and implements the standard OBEX exchange protocol, compatible with PC Windows/Macintosh with OBEX communication programs.

It is possible to disable image transfer with Bluetooth from the control panel (if Bluetooth is not used, it should be disabled to reduce the battery charge and prolong the service life).

7.9 SD Card

There are currently numerous brands and models of SD cards available on the market, not all of them truly meet the requirements. Several recommended SD cards have been tested:

Brand	Sandisk	Kingston	A-Data
Model	SD CARD Sandisk 1GB SD CARD Sandisk 2GB SD CARD Sandisk Ultra II 1GB SD CARD Sandisk Ultra II 2GB SD CARD Sandisk Ultra II 4GB	SD CARD Kingston 1GB SD CARD Kingston 2GB	SD CARD A-Data 4GB class 2 150x

The TakeMS brand SD card should not be used as it does not work with X-POD.

Handheld X-POD recognizes the SD memory card (from 512 Mb to 2 Gb) and SDHC memory card (from 4Gb to 16 Gb).

It is not advisable to use the 4Gb SD card with format SD as it is not standard. If a 4Gb SD card needs to be used, format SDHC is recommended.

The content of the SD Card is made available as a logic unit when connected through USB®. It is strongly recommended not to edit the content of the SD card directly from the PC.

It is strongly recommended to use the iCapture program to transfer files saved on the SD Card. iCapture is software for PC available in the cdrom enclosed with the product. For MACINTOSH platform users, the image synchronization procedure may vary according to the image management software employed.

NOTE: In order to use the MACINTOSH platform, it may be necessary to enable the function to save the images in PNG format, as the program employed may not support Windows WDP format (Hd Photo). This function can be enabled from the control panel. When PNG format is selected, it will take about ten seconds longer to save the images due to the considerable size of the files.

NOTE: It is recommended not to uselessly format the SD card on the PC: the manufacturer's formatting is always the most suitable one. However, it can be reformatted, if necessary, using Windows or Mac; in this case, select FAT or FAT16 (not FAT32) for better performance. FAT32 formatting is required for cards with 4GB or larger capacity.

7.10 USB®

A "micro B" → "type A", highspeed cable should be used to connect handheld _DentalUnitMod2. As soon as the handheld is connected to the PC with the USB cable ® it is forbidden to take x-rays. The user is told that this cannot be done on the handheld screen. Once the USB cable ® is connected it is only possible to synchronize the images.

8 Compatibility with x-ray generators

MyRay sensors work correctly with both conventional x-ray generators called "AC", and newer high frequency generators called "DC".

The following chart specifies the exposure times to be used for two different focal distances (according to the collimator used) between the focal spot and skin; the system considers that the distance between the skin and sensor's active surface is approximately 5 cm.

SETTING EXPOSURE TIMES (IN SECONDS)

	Length of the cone 12" (30cm)		Length of the cone 8" (20cm)	
	0,25	0,16	0,16	0,10
	0,20	0,125	0,125	0,08
	0,16	0,10	0,10	0,063
	0,20	0,125	0,125	0,08
	0,16	0,10	0,10	0,063
	0,20	0,125	0,125	0,08
	0,25	0,16	0,16	0,10

- Tables are referred to a DC generator 60-65 kV, 8 mA, if 4mA is selected it's necessary to double the time.
- If a 70 KV high frequency generator is used, the time given in the table has to be reduced approximately 1/4.
- If edentulous areas are irradiated, the sensor may provide images that are too blackened in the missing areas of the irradiated radiographic subject. In these cases, reduce the time indicated in the table by about 1/4.
- The best results are obtained with a 60 KV high frequency generator, square collimator and 30 cm focus-skin distance (refer to the relevant table).
- To more effectively monitor the distances, the MyRay® centering device equipped with fixed spacer between the centering ring and sensor should be used.
- Before attempting to use the sensor on a patient, practice by taking a few x-ray pictures on unanimated objects with your own x-ray unit.
- Do not exceed the dose in the chart

MyRay® RXDC models HyperSphere or eXTend

This new generation x-ray generator is specifically designed for use with digital systems. In fact, it has 25 different sensitivity settings (“F” factor) to be combined with the digital x-ray receiver used; the sensitivity is selected from the wireless handheld. A 30cm (12”) rectangular collimator, setting F=15, is recommended for use with MyRay sensors. Factor F should be decreased if the images are very dark, according to the positioning technique used.

MyRay® RX AC

This x-ray generator has a selection for digital receivers, shown by a computer icon on the handheld. Make sure the indicator light associated to this function is on. Use the buttons that represent teeth to view the exposure times based on sensitivity factor “S”. This factor can be changed with keys + and - on the handheld: when using with MyRay sensors, set S=8, mode: Sensor

RX DC - CEFLA sc

High frequency x-ray generator prior to the MyRay® RXDC. Even though it has just 3 sensitivity factors for conventional film and a setting for digital receivers, it is still possible to save your favorite exposure times. To correctly use in conjunction with the WDS X-POD system, go to the control panel and make sure 60 kV and DIGIT have been selected. If a long 12” collimator is used, set to 8 mA; on the other hand, with a short 8” collimator, set to 4 mA (the indicator light for the selected function should be on). To save the exposure times suitable for use with the MyRay sensors, select “ADULT patient” and upper molar (18, 17, 16) from the control panel. Use keys + and - to reach the item and set the time to 0.400 s. Lastly, hold down the save key until the confirmation beep is heard.

9 Disposable covers

WARNING!

In order to prevent the transmission of pathogens from one patient to another, it is necessary to always use the disposable covers supplied. Disposable covers are class IIa medical devices and should not be replaced with others having inferior specifications. To obtain more disposable covers, contact the retailer who sold you the sensor. Disposable cover application instructions.

- 1) Place a disposable cover complete with its protective sheet on a flat surface. Introduce the sensor through the opening at one end.
- 2) Push the sensor all the way into the disposable cover, paying attention not to break the transparent material.

- 3) If present, remove the protective support sheet.

- 4) The operation is now completed.
- 5) After use dispose of the covers as “special” waste.

10 Cleaning and disinfecting

Cleaning is the first step of any disinfecting process. Physically scrubbing with detergents and surface-active substances and rinsing with water removes a considerable amount of micro organisms. If a surface is not clean first, the disinfecting process cannot be successful. If a surface cannot be adequately cleaned, it should be protected with barriers.

WARNING!

The handheld external casings and sensor connector are not protected against the penetration of liquids. In any case, disposable infection control sheaths must be used with the sensor to prevent cross-contamination.

Any surfaces that cannot be adequately cleaned should be protected with barriers. The outside of the apparatus must be cleaned and disinfected with a hospital disinfectant labeled as effective against

human immunodeficiency virus (HIV) and hepatitis B virus (HBV) (i.e. low-level) or a hospital tuberculocidal disinfectant o (i.e., intermediate-level) specifically designed for small surfaces.

Follow the manufacturer's operating instructions.

The great variety of medical and chemical products used in a dental practice may damage painted surfaces or plastic materials parts. Tests and studies were conducted to prove that surfaces cannot be entirely protected against etching by all the products available commercially. We therefore recommend to use protective barriers whenever this is possible.

The etching effect of chemicals also depends on surface exposure times.

It is therefore important not to leave the product used on the equipment surface to treat for longer than the working time suggested by the manufacturer.

As the active ingredients found in disinfectants are very harsh, use products that contain maximum:

- **Ethanol at 96%**. Concentration: maximum 30 g. for every 100 g. of disinfectant.
- **Propanol**. Concentration: maximum 20 g. for every 100 g. of disinfectant.
- **Combination of ethanol and propanol**. Concentration: combination of the two should be at most 40 g. for every 100 g. of disinfectant.

The MANUFACTURER has conducted compatibility tests amongst the main disinfectants available on the market and its plastic parts.

The test results indicated that the least aggressive agents are:

- ⚠ Incidin Spezial (Henkel Ecolab).
- ⚠ Omnizid (Omnident)
- ⚠ Plastisept (ALPRO) (not tuberculocide as not an alcohol-based disinfectant).
- ⚠ RelyOn Virkosept (DuPont).
- ⚠ Green and Clean SK (Metasys) (not tuberculocide as not an alcohol-based disinfectant)

The same tests showed that the above-mentioned products can therefore be used by adopting, however, the following preATTENTIONS:

WARNING!



- ⚠ Do not use products containing isopropyl alcohol (2-propanol, iso-propanol).
- ⚠ Do not use products containing sodium hypochlorite (bleach).
- ⚠ Do not use products containing phenols.
- ⚠ Do not spray the selected product directly on the surfaces.
- ⚠ Do not combine the products listed above together or with different liquids.
- ⚠ Any product should be used in accordance with its manufacturer's instructions.

Cleaning and disinfecting instructions

Clean and disinfect with single-use non-abrasive paper (avoid using recycled paper) or sterile gauze.

Sponge cloths or any other re-usable materials are not recommended.

Never spray anything directly on the handheld. Use a dampened soft cloth to clean the handheld.

WARNING!

- In order to clean equipment connected to the power mains, shut off the devices and disconnect the power supply from the outlet before attempting to clean and disinfect the outside.
- Everything used to clean and disinfect must be thrown out after use. Observe current regulations when disposing the material.

Recommendations for waste disposal

Follow the manufacturer’s instructions when disposing whole disinfectant bottles.
Do not let the product into the municipal sewer systems and/or waterways.

11 Material preparation

WARNING!

Fully recharge the battery before attempting to use the X-POD the first time or after it has remained inoperative for a prolonged period.


Depending on the electrical outlets present in the dental facility, select the one that is suitable for the power supply to recharge the battery and insert it as directed in paragraph 14 “Preparing the material”, until you hear a click.

Connect power supply to an adequate electrical outlet. The green light on the power supply will come on. At this point, push the power supply plug straight into the handheld receptacle. If the plug is not already inserted on the cable, put it on pushing it fully into place, aligning the arrow with symbol “-”.

Let the battery charge at least eight hours. The recharging light will light up on the display. A cord is given in the box to more conveniently work. If it is attached to the ring on the back of the sensor the device can be hung around the patient’s neck. To improve convenience, the cord has a snap so that it can be opened in the back.

12 Use of the centering device

The sensor must be held in the correct position using a centering device in order to obtain good x-ray images. A kit containing centering devices specific for front and rear periapical exposures, bite-wing and endodontics is available in the MyRay y® line of products. Use code 97901102 for small sensors and code 97901103 for large sensors to purchase. The individual components in the kit are also available as spare parts. Contact the dealer who supplied the sensor to purchase additional centering devices.

In addition, universal centering devices like RINN® Uni-Grip, KerrHawe® series Bite Sensor or similar can be used. Instructions regarding use of the centering device are given in the Operating instructions included in the centering device kit.

NEVER grasp the sensor with grippers to avoid irreparable damage. Always utilize positioners specifically designed for use with digital x-ray sensors.

Contact the dealer who supplied the sensor to purchase additional centering devices.

Refer to the specific manual supplied by the manufacturer for instructions on how to use the centering device.


13 Possible operating modes

13.1 Stand alone mode

Image transfer through USB connection[®] or WIRELESS BLUETOOTH[®]


Pc with SW iCapture to capture images and iRYS or third party SW for analysis.

13.2 Continuous connection mode

If the handheld's Bluetooth[®] connection is active, any images not yet transferred are saved in the PC's memory whenever it is hooked up to a pc; it is possible to continue capturing images in the mean time.

Image transfer through WIRELESS BLUETOOTH[®] connection

Pc with SW iCapture to capture images and iRYS or third party SW for analysis.

14 Handheld WDS X-POD

14.1 Fast reference guide

- It is strongly recommended to fully charge the battery as soon as you receive the X-POD system. Push the power plug straight into the handheld. Keep the handheld off so that it takes less time to charge the battery.
- The battery is fully charged when the yellow led indicator in the corner stops flashing.
- Once recharged, connect the x-ray sensor with the handheld shut off.
- Make sure the SD card is in the slot.
- Connect the sensor, inserting the connector until you hear a "click".
- Briefly press the power on button to turn on the handheld symbol


- The date and time have to be set when the handheld is turned on the first time. If the handheld does not automatically ask you to do this, remember to set the date and time correctly by tapping

symbol

- Set the desired exposure time, voltage and current for the x-ray unit.
- Position the x-ray sensor using the centering device.
- Start x-ray emission.
- At this point, the x-ray image has been saved on the SD card.
- Briefly press the power on/off button and select  to shut off the handheld.

As far as transferring images to a PC is concerned, consult the paragraphs that deal with the USB®, Bluetooth® and software supplied with iCapture.

14.2 Firmware update from USB®

Handheld X-POD can be updated with software iCapture and the USB® connection. The required procedure is given in the manual that deals with iCapture supplied with the software pack. This can also be done by clicking icon 'info' (

) found in the top right hand corner of the iCapture Monitor card.

14.3 Detailed handheld operation

14.3.1 Setting up the handheld

Once received, the handheld needs to be fully charged before it can be used. See paragraph "Charging and discharging" for instructions on how to charge the battery.

Similarly, the device has to be fully recharged after it has remained inoperative for an extended period (a few months) before it is used again. Recharge the battery with the handheld shut off so that it takes less time to recharge.

Make sure the SD card is inserted (Important information see paragraph "SD card") and that the x-ray sensor is connected (Important information see paragraph "Turning the sensor on and off").

14.3.2 Introduction

Briefly press the power on button to turn on the handheld symbol ).

Whenever the device is turned on, a short diagnostic screen is displayed.

The menu used to set the date and time is displayed when the handheld is turned on the first time.


To set the date and time, press buttons  .

Press  to confirm.

This operation is very important to obtain the exact chronology of the x-rays that will be taken.

The next screen is the main menu from which all the handheld functions can be accessed.

14.3.2.1 Charging and discharging

The device comes with a Lithium-Polymer rechargeable battery. The battery's capacity permits autonomy of a full work day, estimated in over 230 x-ray pictures during 8 hours of consecutive operation. This autonomy is achieved when the battery is fully charged and efficient.

The batteries become less efficient as they get older. It is possible to estimate that its efficiency is reduced to 60% after 500 complete recharge cycles. However, even under this condition, the battery allows about 150 images to be acquired in 5 consecutive hours. When the efficiency of the battery is reduced to an extent at which it is no longer able to sustain the required acquisition rate, call in a qualified technician to change it (original spare part Code 97901336).

While the device is recharging, a symbol similar to this one may appear: .

Vice versa when the battery charge level is low the image of the battery will start fading out  until it starts

flashing. When it begins to flash it means the battery needs to be recharged.

To recharge, plug the power supply provided into an electrical outlet and connect the power supply cable to the handheld. If the handheld is off (recommended), the LED indicator will start flashing until the device has been fully recharged.

If the battery is completely discharged, an alert symbol will be displayed that will not allow the device to be used.

At this point, it is highly advisable to completely recharge the handheld.

It is possible to charge the device with the USB cable[®], however it takes much longer. In this case, the battery receives current from the USB port[®] (5VDC, max 500 mA). In this case, it may take more than ten hours to fully charge. To more quickly charge the battery, use the charger provided and keep the handheld off while charging.

The battery can be charged and discharged hundreds of time however it runs down over time with use. The battery should be changed when the operating and standby times are considerably shorter than usual.

WARNING!

Do not attempt to change the battery on your own. The battery must be changed with a genuine replacement battery code 97901336 by properly trained personnel.

The exact model number of the battery may vary at the MANUFACTURER'S discretion. In this case, the software will be updated by a specially trained technician.

WARNING!

It takes approximately 6 hours to completely recharge the handheld. Use only the FRIWO brand MPP model battery charger provided code FW 7555M/09. Do not use any other battery chargers and do not use battery chargers that are damaged.

The battery charger is supplied with the 4 mm jack already inserted. The jack has a positive pole in the central terminal.

The battery charger mentioned above is available as a replacement part from authorized MyRay dealers.

If the battery is not used for long periods it may still slowly discharge (to keep the date and time). Fully charge after long periods of inactivity before use. In addition, make sure the date and time are correct (set them if necessary).

The handheld is equipped with a timer that automatically shuts it off if accidentally left on for long periods. It will automatically shut off 4 hours after the last image was acquired or the screen was touched.

14.3.2.2 Turning the sensor on and off

In order to acquire images with X-POD connect a MyRay x-ray sensor to the connector found on the handheld.

WARNING!

Plug the sensor connector into the interface making sure to insert it in the right direction. The connector has two small Phillips screws: the screws should face down when the connector is plugged in. When correctly plugged in a click will be heard.


WARNING!

Do not attempt to disconnect the sensor while x-ray pictures are being taken.

To disconnect the connector, press both side tabs to disengage the lock and then remove it without pulling or wiggling to the side.

If the device is turned on, briefly press the power on/off button before unplugging the connector. The operator will have to make a selection on the display, press : the sensor is de-energized, protecting it from an electrical point of view while it is taken out.

The x-ray sensor is expensive and fragile, bothered by static electricity. As a result, always remove it from the interface after the interface has been shut off:

The sensor is also shut off when the handheld is connected to a device via USB® or when the battery is not sufficiently charged. In all these cases, a message is displayed on the handheld advising the user that the unit cannot capture any other images.


If the battery is discharged, the handheld can still work normally and acquire other radiographs if the battery charger is used to supply power.

When the alert symbol given to the side appears, the handheld has detected that the x-ray sensor is not present.

When in this state, all the functions of the device, except for receiving new radiographs, can be used.

If "Sensor Error" appears, shut off the device, unplug the sensor connector, very carefully re-insert it and then turn on the handheld again. Try doing this several times; if the error persists, contact technical service.

14.3.2.3 SD Card

The sensor is also shut off when the SD card is not plugged into the handheld.

The sensor is also shut off when the SD card does not have enough space available to save any other images. In both cases, the alert message shown to the side will be displayed.

 **WARNING!**

Do not attempt to take out the SD card when the corresponding icon  flashes on the handheld display screen. Always keep the SD Card inside the handheld unless it has to be replaced or used to transfer images.

 **WARNING!**

Handheld WDS X-POD is fully operative only when the SD card is inserted. On the contrary, an alert message is displayed. The sensor is also disabled when the SD card is not inserted in the handheld.

 **WARNING!**

If an SD card different from the one supplied needs to be used, it is strongly recommended to format it on the PC. Original formatting by the manufacturer is usually streamlined for that specific SD card.

 **WARNING!**

Insert the SD card with the contacts facing down.

If a 2Gb or 4Gb SD card needs to be formatted, format FAT16 (often indicated only as "FAT") instead of FAT32 is recommended. Access to the memory is faster with FAT16.

 **WARNING!**

Make sure the SD card is not write-protected. Consult the manual that deals with the SD card and turn the selector switch to position "SD not protected".

For safety reasons, the handheld can save and delete x-ray images on an SD-card even if the SD-card is protected as a precautionary measure. It is however good practise to check the position of the protection mark and leave it as shown in the figure. The user has to confirm before any images are actually deleted.



14.3.3 Capturing images

Once the previous operations have been completed, it is possible to start capturing images. Handheld X-POD is designed to save a number of tests or images that belong to the same patient or imaging session. Each test (or series) is placed in a separate folder. Each folder can hold up to 25 images. Each test (or series) is automatically frozen when a new one is opened. Once a test (or series) is closed, it is no


longer possible to add any images to it. First of all, tap icon  and enter the name of the test (or series of x-rays to be taken). **Before attempting to take an x-ray picture, always verify the device's LED is BLUE and flashes**. As soon as the sensor receives the x-rays, the image is saved with its own name on the SD Card.

The captured images can be consulted later on by pressing  and selecting the desired name (paragraph "Managing and consulting the images").

Tap icon  to access the tables with typical exposures. Tap the figures to scroll the pages.

Tap icon  to go to the Settings menu (paragraph "Settings menu").

As far as transferring images to a PC is concerned, consult the paragraphs that deal with the USB®, Bluetooth® and software supplied with iCapture.

Put the sensor into place using the centering devices provided. Compatible x-ray generators should be used. Always follow the precautions and recommendations given in the dedicated chapter.


X-POD is not able to capture new x-ray images while it is connected to a PC or other device via USB®; in this case, do not emit any x-rays.

The acquired images are automatically saved on the SD card, in the native format not included. The number of images that can be acquired depends on the capacity of the SD card used

If the SD memory card is full, copy the images in a PC or other storage unit before deleting them or replace the SD card with a new one.

If a "Generic Error" is signaled with the symbol shown to the side during any operation, try to reset the device and perform the desired operations again. If the problem persists, consult paragraph "Troubleshooting". If the "Generic Error" is still present even after consulting the "Troubleshooting" paragraph, contact technical service.

14.3.4 Managing and consulting the images

Tap on the main screen to consult the series of captured images (image to the side). Tap to scroll the entire list of series saved on the SD CARD, if there are more than seven. It is also possible

to find a certain name by selecting and entering the desired series.

Tap the desired name twice to display the images of a series. The series of images already closed are blocked by default and other images cannot be added. They can be released. To do this, press key

'lock closed' which will become 'lock open' .


When the lock is opened, the name of the series currently selected can be changed and images can be added to the series. The selected series will be highlighted in red, to remind the user that a previously closed series has been opened. The lock will close and the series will be blocked again whenever an unlocked series is quit.

To change the name of a series, press key 'rename' and enter the new name. After entering the new name, press "OK" to confirm. Once all the images of a certain series have been transferred to the

PC, the 'rename' icon will become grey , as the command can no longer be used.

It is always possible to delete a series of images provided that the user confirms they are to be canceled. To delete a series, tap icon . In the event that a series was not completely transferred to the PC the warning message will have a note to remind the user that all the images have not been saved in the PC.


If a few images are present, a mosaic display of the preview will appear. To select a single image, simply tap it. As a result, a fully screen image will be displayed.

Functions “Pan” and “Zoom” () are also integrated. There are three zoom levels: touch the area on the screen to be enlarged to obtain a full screen image. After tapping three consecutive times a full screen image is displayed. To drag the image, keep your finger on the screen and move the image in the desired direction.

If there are several images in the series, scroll by pressing them one after the other:

When a corner of an image is light blue, it means that the image has already been transferred to the PC.

If an image is not present in the selected position, a screen with a blue cross indicating the image is missing will be displayed.

If an x-ray picture is acquired at this time, it will be saved in this position, thereby replacing the missing image symbol.

This operation can be repeated several times, if necessary.

14.3.4.1 Editing and measuring images

Whenever an image is selected the icons are always arranged in an intelligent manner regardless of how the handheld is tilted.

Once an image is displayed, a number of tasks can be performed and adjustments made to it.

Activates and deactivates the menu for adjusting the images. Zoom and pan are active in this mode.

Adjusts brightness.

Tap symbol or on the touch screen respectively to increase or decrease brightness.

Adjusts contrast.

Tap symbol or on the touch screen respectively to increase or decrease contrast.

Adjusts the range.



Tap symbol or on the touch screen respectively to increase or decrease the range.

Image measurements.

Measurements can be taken on the x-ray image. If key "measurement points" (highlighted in green in figure) is active, measurement points can be entered. The first time the screen is tapped indicates the point at which the measurement starts; each subsequent tap will be summed to the length of the section thereby obtained. If key "measurement points" is pressed again, it will

be deselected () and the position of the individual points can be modified. The individual points can be moved by dragging the round targets.

It is possible to calibrate the measurements on the images by tapping icon . Once tapped, this icon will open the menu in which the length of the known measure can be entered. Once calibration has been set, the gauge symbol will appear to the left of the measurement. Calibration can also be eliminated.

Press key to cancel the measurement points on the screen. All the points will automatically be canceled.

NOTE: In order to quit the measurement screen, deselect measurement mode by opening the tool menu and tapping the ruler. As a result and only in this case it is possible to go back to the previous menu used to manage the series of images.


Undo changes.

Click this icon to undo all the changes made. The image will go back to its original appearance.

Image/teeth correspondence

NOTE: Dental arch ("Dentition Chart") activated. This screen is normally accessible; if you want to disable it, see page 70.

Select this icon to view the upper/lower arches (use the icons shown in the figure to select). Select the teeth from the screen to associate them to the x-ray. The selected teeth will be highlighted in yellow. The letter 'R' indicates the patient's right side. The main properties of the displayed image are highlighted on the bottom. The user can access the menu used to orientate the images from this screen by pressing icon .

The user can access the menu used to orientate the images from this screen by pressing icon .

Deletes selected image.

Press this icon to delete an image saved on the SD CARD. The user has to confirm the image is to be deleted.

14.3.4.2 Modifying image orientation

Once an image has been selected, the orientation with which it is to be displayed once downloaded in a PC can be changed.

Tap the icon on the top left-hand side to access the orientation menu.



The yellow arrow on the tooth icon  identifies the “top” of an image.

If the tooth is yellow, it means that it has been assigned an “image/ tooth” correspondence.

When the handheld is turned, a light blue arrow will blink for a few seconds to suggest the new orientation. Simply tap the icon to accept the suggested one.

If the icon is tapped when at rest, i.e. not blinking, it is possible to manually set the desired orientation. This feature is convenient for professionals who work with the handheld placed on a horizontal surface. The screen to the side will appear. Once orientation has been manually set, tap the check mark to quit the screen and save the new orientation.

Tap the icon in the middle to quickly reach the screen used to select single teeth.

14.3.5 Settings menu

Tap icon  to go to the Settings menu. The screen shown to the right with the possible settings will be displayed.

 Setting the date and time.

Press buttons  and  to set the date and time.

Press  to confirm.

 Touch screen setting.

If the touch screen is not well set, a simple procedure can be carried out to set it. Tap the icon to set up and then follow the indications of the arrow and touch the given spot.

 Bluetooth settings



The name of the handheld is given at the top of this screen. It can be changed by clicking it and then entering the new name.

Confirm with "OK".

If the Bluetooth connection requests the PIN it can be entered on the second line.

It is possible to enable or disable the Bluetooth connection, as desired, by pressing "ON" or "OFF" respectively.

In addition, it is possible to set compression for image transfer; touch the key under the hare (illuminated dots) to speed up image transfer to the detriment of resolution.

Format and filter setting.

On this screen you can activate or deactivate saving the images in PNG and JPG format.


To activate the desired function tap the ON icon, and to deactivate it tap the OFF icon.

Sharpness

This setting adjusts the filtering applied to the image during acquisition; it therefore has no effect on previously acquired images.


The native images in PNG and WDP format are not modified.

As the images in PNG and WDP format are native and hence not processed, they are not altered with this filter; they are however marked with the filtering information used so as to inform any image management programs installed on the PC.

Generic settings.

On this screen you can activate or deactivate the sounds (BEEP), the rotation function of the images and the dental arch ("Dentition Chart").

Information about X-POD:

Screen with system information.

The first line at the top shows the software version installed in the handheld ("FW VER: XX.XX").

Consumption.

It is possible to set three different energy savings modes. When all three batteries are on, the handheld will use more energy; the display will remain fully illuminated longer each time it is touched.

If fewer batteries are on, the handheld's autonomy will be much longer but the display's brightness will be reduced very quickly.

14.3.6 Indicator LED

A multicolor led indicator is present in the corner of the handheld. It shows when the sensor is ready to receive images. **The only thing really important to remember is that the handheld is “ready” to receive x-rays when the led indicator is blue (steady, or flashing slowly if in standby).**

Color of led	Handheld status.
Off	The handheld is fully de-energized.
Flashing blue light	Handheld operative, ready to receive images.
Flashing yellow light, with handheld on	The handheld is not ready to receive images. The sensor is either not correctly installed or is being calibrated.
Flashing yellow light, with handheld off	The handheld is not ready to receive images. The battery is being charged with the external charger or USB port [®] . The light will stop flashing after the battery has been charged.
Red flashing light	Hardware failure. Contact technical service.

14.3.7 Image transfer

The images acquired with the handheld can be transferred to another device (PC, printer, etc..) in 3 ways:

- USB[®] connection
- Bluetooth[®]
- SD card

USB[®]

The USB cable connector[®] must be plugged into the handheld connector.

Use the USB cable provided or a similar cable that is not more than 2 meters long.

Make sure the cable is marked “USB 2.0 high speed”. If the quality of the cable is not adequate, image transfer will take much longer

A Micro-B (not Mini) connector is used, increasingly employed for next generation handhelds and mobile phones. In this case, the images are very quickly transferred.

The transfer mode depends on the device to which the images are transferred.

For example, if the dedicated software is used, transfer to the PC takes place as soon as the two devices have been synchronized whereas if no software is used, the handheld appears as an external removable disk similar to the USB[®] Flash memory pen drive.

When handheld WDS X-POD is connected to the pc or any other device via USB[®] the sensor is always off (see figure to the side). It is automatically turned back on when the device is disconnected.

Bluetooth[®]

If the dedicated software is used on the PC and the Bluetooth[®] connection is on, the content of the SD card is automatically synchronized with the PC hard disk. This can also be done simultaneously with image acquisition.

The Bluetooth[®] connection is active when the relative symbol is white: . If the symbol of the Bluetooth[®] connection flashes, it means that an image is being transferred to the PC.

The Bluetooth connection[®] is NOT active when the relative symbol is grey: .

Loss of the Bluetooth signal[®] does not affect acquisition of the x-rays taken simultaneously. The acquired images can be transferred when connection is restored or in one of the other modes. A Bluetooth dongle is provided to be used in the PC for Bluetooth connection (consult the software manual for installation instructions).

SD card

The images are saved by handheld WDS X-POD on a SD card which is put in and taken out of the slot.

The content of the SD card can be viewed and transferred with a standard SD Card reader (not supplied).

14.3.8 Handheld display screen

The handheld has a color "Touch Screen".
The captured x-ray images can be retrieved from the graphic interface

menu by pressing  and then selecting the name of the desired series of x-rays pictures.

After tapping the name, the preview of the captured images is shown directly on the display (see figure to the side).

The images displayed for previewing cannot be used for diagnostic purposes but to simply verify if the sensor is in the correct position in relation to the area to be examined.

In order to more accurately examine these images, they should be displayed on a large screen and interpreted by a dental professional.

 If the display breaks, the touch-screen membrane holds back any fragments of glass and the handheld is no longer useable. Promptly return to the service center for repair.

15 Quality of the x-ray images

Unlike what usually occurs when regular radiographic film is used, systems that adopt MyRay sensors tend to automatically correct any exposure errors, providing images that are in any case usable.

Notwithstanding the dynamics of the sensor allow images with a wide range of grey levels to be acquired, common PC monitors show only 256. As a result, in most cases the software will be able to extract a good quality image even if exposure is not correct. However, there are limits beyond which good results cannot be recovered.

The x-ray sensor is more sensitive than x-ray film therefore exposure times usually have to be reduced.

See the relevant instructions given in paragraph **Compatibility with x-ray generators**.

To obtain top performance of digital x-ray sensors, it is important to keep in mind that there are some differences compared to film. Insufficient exposure is clearly seen on x-ray film as the areas corresponding to soft tissues are dark. On the other hand, when a digital sensor is used, the image's background noise increases (salt and pepper noise) and the tonal range is insufficient.

Excessive exposure (time too long) on x-ray film causes the image to be too dense (dark) while the image will lose contrast in a digital sensor.

It is a common mistake to confuse excessive exposure with insufficient exposure therefore further increase the exposure times.

It is important to check and note this limit with your own x-ray system in order to be sure not to exceed it during dental treatment as the quality of the images obtained under these conditions will be poor or even unusable.

 **WARNING!**

Before attempting to take x-ray pictures on patients, it is advisable to take some test pictures on phantoms, comparing the results obtained to the usual ones. Identify the best exposure conditions for your own x-ray system through trial and error.

16 Problem solving

Problem found	Possible causes	Remedies
Doubts concerning the sensor efficiency.	Falling, banging, malfunctioning.	Do not use the sensor on a patient. Conduct tests by capturing x-ray images using a phantom.
Image lost during transfer with Bluetooth® connection.	Management software or PC error, input power cut.	It is possible to recover acquired images directly from the SD card or through the USB connection®.
An ERROR message followed by a number appears on the PC.	Sensor malfunctioning.	Do not use the sensor, contact the Technical Service Center.
USB connection® does not work properly.	Defective cable, PC USB® port does not work.	Use the USB HUB® provided and repeat the correct operations. Try using different USB ports® in the PC. If the problem persists, contact technical service.
An image is acquired but with poor tone range and/or massive background noise.	Underexposed image.	Use a longer exposure time, make sure the x-ray generator works correctly.
The image is acquired, but the dark parts invade other areas, making the image unusable.	Over-exposed image.	Use a shorter exposure time, check the settings of the x-ray generator.

A. FCC Regulatory information

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device contains FCC ID: POOWML-C40.

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